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
Day Two: Use a rotary evaporator to remove the solvent. Using a 5.75" Pasteur pipette, prepare a microscale column with neutral alumina as the absorbent and dichloromethane as the eluent. (See pages 678-685 of your laboratory textbook for background.) Place a small piece of cotton in the pipette and carefully push it to the bottom. The cotton should allow liquid to freely move through the column. Fill the pipette with approximately 1.5" of alumina. Add dichloromethane (2 mL) to the column allowing it to drain through the alumina until the solvent surface is just above the alumina surface. This fraction of dichloromethane can be properly discarded. Prepare a 25-mL round-bottomed flask to collect the product. Transfer the product from the 250-mL round-bottomed flask to the wet column. Rinse the round-bottomed flask with dichloromethane (1-2 mL) to dissolve any remaining product and add this to the column. Allow the mixture to flow through the column until the solvent level just reaches the alumina. Add dichloromethane

(2 mL) to the column and again allow the mixture to pass through the column as before. Repeat with another portion of the solvent (2 mL) to elute the product from the column. Use a rotary evaporator to remove the solvent.

Day three: Characterize the starting material/product mixture by  $^1\text{H}$  NMR. Estimate the amount of starting material that is present. (You may want to get a NMR spectrum of the starting material for ease in comparison.) Take an optical rotation and calculate the percent enantiomeric excess based on the concentration of the alcohol present in the sample.

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